

# Rascal P327



**Mid-wheel powerchairs**

**Your new route to independence**

# Rascal P327

## Rascal P327 Specification:

Type Class	B / 2
Overall length	112cm / 44"
Overall width	69cm / 27"
Floor to seat height	54- 56.5cm (min-max)
Seat width	47cm / 18.5"
Seat depth	40.6cm / 16"
Maximum carrying capacity	136kgs / 21 stone
Battery voltage and capacity	12V / 55Ah x 2
Maximum speed	6.4km/h / 4mph
Ground clearance	5cm / 2"
Maximum kerb height	5cm / 2"
Range*	up to 32km / 20 miles
The maximum safe slope	6°
Turn-around width	1.20m / 4'
Colour	Red / Teal

\*Range on full charge and flat ground.

## Main Features:

- Flip-up fully adjustable armrests
- Full suspension
- Swing-away joystick
- Convenient joystick charging port
- Adjustable captains seat
- Flip-up footplate



## Rascal P327 Profile:

The Rascal P327 is a 4 mph Mid-Wheel Drive outdoor/indoor powerchair, with durable full suspension and comfortable cushioned Captain's seat. The P327 delivers great performance with its 55Ah batteries contributing to its excellent range, outstanding manoeuvrability and a smooth ride. Enquire today about upgrading to a larger seat at 20.5". P327 Mini, P327 XL and Seat Lift models also available.

### Full Suspension



Up to 32km




Local Rascal Mobility Retailer:

Max. Carrying Capacity  
136kg / 21 Stone



Maximum Speed



Specifications subject to change without prior notice  
E. and O.E- v 4 Final-July 2018

