

**Pride®**

**FINESSE SPORT**



Introducing the compact 8 Mph Apex Finesse Sport scooter from Pride Mobility. Combining a high-torque motor package and Comfort-Trac front and rear suspension the Finesse Sport is the perfect scooter for the most active individual.

With hi-visibility front and rear LED lightning, exceptional leg room, easy to use delta tiller controls and upgraded 50 Ah batteries the Finesse Sport is the ideal scooter for any journey.



@Pridemobilityuk



Version 1.0 | September 2018



32 Wedgwood Road, Bicester, Oxon, OX26 4UL | Tel: 01869 324600 | Fax: 01869 323070 |  
Web: pride-mobility.co.uk | Email: sales@pride-mobility.co.uk



# Finesse Sport | Standard Features



Front and Rear Suspension



Front and Rear LED Lighting



Easy to Use Delta Tiller Controls



Cabin Light



Alloy Style 10" Wheels



Comfortable Hi-back Seat

# Finesse Sport | Colour Options














Candy Apple Red



Viper Blue



## Technical Specifications

Weight Capacity	Max Speed	Ground Clearance	Turning Radius	Overall Length	Overall Width	Front/Rear Tyres	Total Weight w/o batteries	Heaviest piece	Standard Seating	Batteries Required
 28 St 8lbs. (181 kg)	 8 Mph (12.8 km/h)	 5.6 cm <sup>4</sup> (2.2 inch's)	 156.6 cm (61.25 inch's)	 119 cm (47 inch's)	 64.8 cm (25.5 inch's)	 Front: 25.4 cm Rear: 25.4 cm	 11.5 Stone (73 k.g.)	 5.5 Stone (35 k.g.)	 Hi-Back Seat 13 k.g	 x2 50 Ah deep cycle

### Warranty

1. Life time warranty on frame
2. Two-year limited warranty on drive train including: transaxle, motor and brake.
3. One year limited warranty on all electrical parts
4. 12 month warranty for batteries

### Notes

1. Speed varies with user weight, terrain type, battery charge, battery condition and tyre pressure.
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.
3. Battery weight may vary based on manufacturer.
4. Measured at motor
5. Information correct at time of printing
6. Images used for illustration purposes only.